



Expertise In Research Anesthesia, Equipment And Supplies

## **Imaging Chamber 78919969**

Patterson's five-mouse Imaging Chamber is designed for use with the Caliper/Perkin Elmer IVIS Optical Imaging System. The Imaging Chamber is sealed such that WAG's do not escape into the workplace, and, with the use of Hepa Filters, pathogens cannot get in or out of the Imaging Chamber. The five-mouse manifold increases throughput for studies using the Optical Imaging System.

**78920341** – Imaging Chamber with three-mouse Ivis manifold

