

EquaFlow Manifold



The leader in anesthesia
equipment and supplies.

Manifolds designed to reduce exposure to waste anesthetic gas

Our EquaFlow manifolds are designed to prevent unwanted animal awakenings, scavenge 99.99% of waste gas and enhance animal ergonomics. Available in three sizes, these manifolds can be used in the prep area outside an optical imaging station, in our 3 or 5 mouse-sealed imaging chambers for Bruker or Perkin Elmer or inside an imaging system as a stand-alone manifold. The EquaFlow manifolds require attenuated in-house vacuum or Patterson's Evac 2/4 to operate.

EquaFlow 3I
78919271

EquaFlow 3XL
78919266

EquaFlow 5
78919268

EquaFlow Low Profile
78926554



EquaFlow 3I



EquaFlow 3XL



EquaFlow 5



EquaFlow
Low Profile

We can help you build your system.

Call our experts today to place your order: 800.877.8989